



GE Vivid 3 - Fully Refurbished

The GE Vivid 3 ultrasound system is perhaps the most versatile piece of ultrasound equipment you can buy. With high speed and high quality imaging, the GE Vivid 3 ultrasound can handle up to 18 probes, making it perfect for almost any application. And despite its flexibility and power, it's still amazingly easy to use. The GE Vivid 3 brings ultrasounds to a new level. This machine, like many other GE ultrasound machines, uses the TruScan Architecture imaging platform. TruScan provides physicians with excellent image manipulation and storage, not to mention its superior ease of use. The Vivid 3 will improve the radiology department with its excellence in image quality, productivity, and patient care.

Features:

- DICOM connectivity
- Automatic Tissue Optimization (ATO)
- 18 transducers for a wide range of applications
- Advanced 3D option for vascular studies
- Enhanced OR performance with TEE and intra-operative probes
- Stress Echo options for faster exams that allow for protocol changes mid-exam
- Includes C721 transducer for pediatric and neonatal applications
- Coded Harmonics Image Quality
- TruAccess Total Data Management
 - MPEGvue option allows fast saving of patient records onto CDs
 - Utilize custom report templates by the choosing from a variety of formats and layouts
- Automatic Tissue Optimization (ATO) makes it easy to get quality images quickly
- Anatomical M-Mode enables off-axis correction when the heart isn't positioned as expected or is not the expected shape
- Dimensions: 52" (H) x 25" (W) x 44" (D)
- Weight: 350 lbs
- 15" CRT color monitor

